

Breathpal Reproducibility over Time and after Repeated Connection

© 2023 Ben Bulsink

Version 1.0-March 17, 2023

Content

Background.....	1
RIP Measuring principle.....	1
Test setup.....	2
Measurements	2
Discussion.....	3
Long term stability.....	3
Disconnection and reconnection	3
Changing of straps.....	3
Conclusions.....	4

Background

The Breathpal RIP sensor can be used to measure the breathing movements of the body.

To use the sensor over a longer time and between separated sessions, it is important to know the stability of the reading over time and between separate attachments. This document describes the method of investigation of this stability and the results.

RIP Measuring principle

The RIP measuring method actually measures the surface (cm²) within the circumference of the strap that is spanning the body. This surface is related with the tidal volume (cm³) of the breathing, with a beforehand unknown factor, and with various unknown stable and less stable errors. See <https://breathpal.nl/wp-content/uploads/2021/12/Breathpal-measurement-quality-report-vs-1.4.pdf> for an overview of errors in the translation from tidal volume to surface.

Test setup

The tests are executed by using the Plastic Lung, a device supplied by MST (Figure 1). It allows to put straps around two halfpipes, and precisely vary the distance between the halfpipes, and so change the surface within the strap loop with a reproducibility (scale division) of 0,16 cm².

The sensitivity of the sensor is set to be 4.6 Arbitrary Units (AU) per cm² surface change (Plastic lung reading: 17,5 cm., enclosed surface 472,5 cm²)

Measurements

The sensor is set up, using the plastic lung, as in below picture. The distance between the half pipes is locked to a fixed value of 17.0 cm.



Figure 1: Setup of long term stability measurement

Between each measurement, the USB powering and the abdominal strap had to be disconnected while the battery runs empty over one day, which causes that repowering does not initialize USB connecting well then.

Below the (AU) readings over a time span of three weeks.

Four strap are used:

A: 74 cm length, in rest

B: 76 cm length, in rest

C: 72 cm length, in rest

D: 68 cm length, in rest

Table 1: Readings over longer time

	Thorax reading strap A (AU)	Abdominal reading strap B (AU)
16-2-2023 11:48	7418	2958
17-2-2023 13:00	7417	2960
09-03-2023 10:40	7420	2961

Now, on date 17-03-2023, repeatedly both straps are disconnected and reconnected, sometimes with some manipulation. The remarks and readings are listed in below table.

Table 2: readings after reconnections

Measurement number and remarks	Thorax reading strap A (AU)	Abdominal reading s trap B (AU)
1	7428	2970
2	7428	2970
3a	7426	2969
3b after some shifting and stretching	7428	2967
4	7429	2966

Now, on 17-03-2023, the readings are taken, and the straps are exchanged and changed by two others.

Table 3: Readings with different straps

Measurement number and remarks	Thorax reading (AU)	Abdominal reading (AU)
After power-up and connecting	7411 strap A	2956 strap B
After powerdown and reconnecting	7412 strap A	2955 strap B
Exchange straps	7440 strap B	2926 strap A
Use other straps	7409 strap C	2676 strap D
And exchange these	7151 strap D	2921 strap C

Discussion

Long term stability

The long term stability of Breathpal is very good. For making readings over a longer time, it was decided to switch the computer off, which made it necessary to disconnect and the straps before each reading (in order to avoid improper USB enumeration, due to technical choices in Breathpal). With as little manipulation of the textile straps as possible, the readings had a maximum deviation over 3 weeks of only 3 AU (0.65 cm²). A better procedure would be to keep the reading computer up and running over the long time, which would give a more precise and reliable figure of the long term stability. However, in practical use, this situation is unlikely to appear.

Disconnection and reconnection

The disconnection and reconnection of the Breathpal straps also causes different readings on the same mechanical situation. Actually, the readings in Table 1 also needed disconnection and reconnection of the straps. Table 2 lists these effects where more manipulation is done. Some mechanical shifting and stretching of the straps does not have much influence.

Changing of straps

The effect of changing straps by other versions (same material, but slightly different lengths, is listed in Table 3. Changes of 250 AU are found, just by a few cm strap rest length difference. If absolute value readings need to be compared over some sessions, it is necessary to use the same straps for these sessions.

Conclusions

The Breathpal electronics show very good stability over time. The strap material is based on soft tissue, which can easily change the surface within the copper wire connector which it supports. Therefore slight manipulations of the straps can change the sensor reading. The readings do largely change between straps that have different length in rest.